

TECHNICAL DATA SHEET: HYDROLYZED KERATIN PROTEIN

Page 1 of 2

DATE OF PREPARATION: May - 4<sup>th</sup>-2015  
VALIDITY: May - 4<sup>th</sup>-2016

VERSION 02

BATCH DELIVERY DATE (dd-mm-yyyy)	LOT No./SAMPLE No.
20-04-2015	002/001

PRODUCT	
NAME	LIQUID HYDROLIZED KERATIN PROTEIN

CHEMICAL ANALYSIS	
TEST NAME	Results
Organoleptic evaluation	Amber colored liquid
Residue on ignition (ash)	≤ 2 %
Protein	≥20 %
pH (Solution 5%)	6.0-7.0
As	≤ 0.5 ppm
Cd	≤ 10 ppm
Cr	≤ 10 ppm
Hg	≤ 0.5 ppm
Ni	≤ 10 ppm
Pb	≤ 1.0 ppm
Total enterobacteriaceae	≤ 1000 CFU/g
Yeasts and molds	≤ 100 CFU/g
Escherichia coli	≤ 30 MPN/100g
Salmonella	Absent
Staphylococcus aureus	Absent
Origination	Hoofs and horns

TECHNICAL DATA SHEET: HYDROLYZED KERATIN PROTEIN

Page 2 of 2

**Handling and storage:**

Keep container tightly closed, storage should be 2-25°C and dry. Just for analytical purposes.

<b>MADE BY</b> ORBI analyst	<b>APPROVED BY</b> ORBI Production Plant - Chief
Victor Villa	Andrés Mejía
	